# Super Sani-Cloth®

GERMICIDAL DISPOSABLE WIPE  
LINGETTE GERMICIDE JETABLE

EXTRA LARGE WIPE  
TRES GRANDE  
LINGETTE

Effective against / Efficace contre :

**BACTERIA**

- *Bordetella pertussis* • *Burkholderia cepacia* • *Campylobacter jejuni* • *Escherichia coli* [E.coli] • *Escherichia coli* O157:H7
- *Klebsiella pneumoniae* • *Pseudomonas aeruginosa* • *Salmonella enterica* • *Staphylococcus aureus*

**MULTI-DRUG RESISTANT BACTERIA**

- *Acinetobacter baumannii* • *Enterobacter cloacae* NDM-1 positive • ESBL Producing *Escherichia coli* [E. coli] • *Klebsiella pneumoniae* - KPC-2 positive, ST258 • Methicillin Resistant *Staphylococcus aureus* [MRSA] • Vancomycin Resistant *Enterococcus faecalis* [VRE]

**VIRUSES\***

- Adenovirus Type 5 • Herpes Simplex Virus Type 2 • Human Coronavirus Strain 229E • Influenza A/Hong Kong • Influenza A [H1N1] virus • Respiratory Syncytial Virus [RSV] • Rhinovirus • Rotavirus Strain WA • Vaccinia virus • Kiltis Pandemic 2009 H1N1 Influenza A virus

- BLOODBORNE PATHOGENS** • Hepatitis B virus [HBV] - Duck HBV\* • Hepatitis C virus [Human] [HCV] - Bovine Diarrhea Virus\* • HIV-1 [AIDS virus]

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- TB** • *Mycobacterium bovis* - BCG [tuberculose]

- CHAMPIGNONS PATHOGENES** • *Candida albicans*

BACTERICIDAL, TUBERCULOCIDAL, AND VIRUCIDAL\* IN 2 MINUTES  
BACTÉRICIDE, TUBERCULOCIDE, VIRUCIDE\* EN 2 MINUTES

**ACTIVE INGREDIENTS:**

n-Alkyl (68% C <sub>12</sub> -32% C <sub>14</sub> ) dimethyl ethylbenzyl ammonium chlorides.....	0.25%
n-Alkyl (60% C <sub>12</sub> -30% C <sub>14</sub> , 5% C <sub>16</sub> , 5% C <sub>18</sub> ) dimethyl benzyl ammonium chlorides.....	0.25%
Isopropyl Alcohol.....	55.00%
<b>OTHER INGREDIENTS</b> .....	44.50%
<b>TOTAL</b> .....	100.00%

Does not include the weight of the wipe.

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chlorures de n-alkyl (68% C <sub>12</sub> -32% C <sub>14</sub> ) alkyldiméthyl (éthylbenzyl) ammonium.....	0,25%
chlorures de n-alkyl (60% C <sub>12</sub> -30% C <sub>14</sub> , 5% C <sub>16</sub> , 5% C <sub>18</sub> ) alkyldiméthyl (éthylbenzyl) ammonium.....	0,25%
Alcool isopropyle.....	55,00%
<b>AUTRES INGRÉDIENTS</b> .....	44,50%
<b>TOTAL</b> .....	100,00%

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8K086900

0100310819000581

\*KILLS HIV-1 (AIDS VIRUS), HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) ON THE PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS AND PRE-CLEANED EXTERNAL SURFACES OF ULTRASONIC TRANSDUCERS AND PROBES PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in two (2) minutes at room temperature (68-77°F) in healthcare or other settings in which there is an updated likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surface/objects likely to be soiled with blood or body fluids can be associated with the potential transmission of Human Immunodeficiency Virus Type 1 (HIV-1) associated with AIDS, Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV).

**Areas of Use: Hospital and Healthcare Settings:** Ambulatory Surgical Centers (ASC), Clinics, Dental Offices, Dialysis Clinics, Home Health Care, Hospices, Hospitals, Laboratory, Nursing homes, Physical Therapy, Physician's offices, Radiology, Rehabilitation, Transport vehicles, Critical Care Areas: CCU, Emergency Rooms, ICU, Neonatal Intensive Care Units (NICU), Operating Rooms, Pediatric Intensive Care Units (PICU), Surgery and Surgical Intensive Care Units (SICU), Hospital, Healthcare, and Critical Care Use Sites: May be used on hard nonporous surfaces of bed railings, blood glucose meters, cabinets, carts, chairs, counters, dental unit instrument trays, exam tables, gurneys, iceboxes, IV poles, stethoscopes, stretchers, tables, telephones, toilet seats and hard nonporous outside surfaces of ambulances and dental curing lights; diagnostic equipment, patient monitoring equipment, patient support and delivery equipment. This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

**DIRECTIONS FOR USE:** It is a violation of Federal law to use this product in a manner inconsistent with its labeling. To Dispense Wipes: Remove lid and discard inner seal from canister. First corner of wipe roll, remove first wipe for use, best corner of next wipe into a point and thread through the hole in the canister lid. Pull through about one inch. Replace lid. Dispense remaining wipes as necessary by pulling out at an angle. When not in use keep center cap of lid closed to prevent moisture loss. **TO DISINFECT AND DEODORIZE:** To disinfect nonporous contact surfaces only. Unfold a clean wipe and thoroughly wet surface. Allow treated surface to remain wet for a full two (2) minutes. Let air dry. For heavily soiled areas, use a wipe to pre-clean prior to disinfecting. These directions also apply to *Mycobacterium bovis* BCG (Tuberculosis) at 68°F (20°C).

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.**

**Personal protection:** When using this product, wear disposable protective gloves, protective gowns, masks, and eye coverings when handling HIV-1 (AIDS Virus), HBV or HCV infected blood or body fluids. **Cleaning procedure:** All blood and other body fluids must be thoroughly cleaned from surfaces and objects before disinfection by the germicidal wipe. Open, unfold and use first germicidal wipe to remove heavy soil. Dispose of infectious materials: Used wipe, blood and other body fluids should be disposed of according to local regulation for infectious waste disposal. Contact time: Use second germicidal wipe to thoroughly wet surface. Allow to remain wet two (2) minutes, let air dry.

**PRECAUTIONARY STATEMENTS** Hazards to Humans and Domestic Animals

**WARNING:** Causes substantial but temporary eye damage. Do not get in eyes or on clothing. Avoid contact with skin. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using medicine. Remove and wash contaminated clothing before reuse. **FIRST AID:** Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. If on skin: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. **PHYSICAL OR CHEMICAL HAZARD:** Combustible. Do not use or store near heat or open flame. Do not use on natural marble, windows, unpainted wood, brass, clear plastic or colored glass. Test wipe on small inconspicuous area first. **STORAGE AND DISPOSAL:** Do not contaminate water, food, or feed by storage and disposal. Storage: Do not store near heat or open flame. When not in use keep center cap of lid closed to prevent moisture loss. **Toilette Disposal:** Do not reuse. Toilette: Dispose of used toilette in trash. Do not flush in toilet. **Container Disposal:** Nonrefillable container. Do not reuse or refill this container. Offer for recycling, if recycling is not available, put in trash collection.

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