



YOU ARE THE FRONTLINE OF INFECTION PREVENTION.




INTERVENTIONAL CARE PATIENT CARE ENVIRONMENT OF CARE

ES HAS A POWERFUL PARTNER IN INFECTION PREVENTION.

Sani-Cloth® products are specifically designed to support Environmental Service Professionals in their vital role to help stop the spread of preventable infections, whatever the challenge.

Features

- Pre-moistened, ready to use
- No measuring, no mixing, or waiting
- Color-coded for ease of implementation, replacement and compliance
- Compliant with CDC Guideline Disinfection and Sterilization in Healthcare Facilities and AHE Practice Guidance for Environmental Cleaning¹
- Designed for Environmental Services convenience & efficiency with bilingual materials



#1
DISINFECTANT WIPE IN HEALTHCARE²

¹Practice Guidance for Environmental Cleaning, Association for the Healthcare Environment, 2012. Guideline for Disinfection and Sterilization in Health Care Facilities, Centers for Disease Control and Prevention Healthcare Infection Control Practices Advisory Committee, 2008.

²GHX 2018 Q3 report.



Sani-Cloth® Bleach Germicidal Disposable Wipe

Ideal for disinfecting high risk areas endemic with *Clostridium difficile* spores and Norovirus.



Sani-Cloth® AF3 Germicidal Disposable Wipe

Alcohol and fragrance free formula is ideal for use around patients, residents and staff, especially those with respiratory sensitivities.



Super Sani-Cloth® Germicidal Disposable Wipe

Ideal for daily use in fast-paced environments that require short contact times and broad coverage of microorganisms.





Sani-Cloth® Prime Germicidal Disposable Wipe & Sani-Prime® Germicidal Spray

Powered by ACCELOQUAT™ formulation, a next generation proprietary blend of quaternary ammonium, isopropyl alcohol (IPA), and ethanol.



Easy Screen® Cleaning Wipe

A compatible touchscreen cleaner that works at the speed of technology, featuring the power of 70% isopropyl alcohol (IPA).



Sani-Cloth® Point-of-Care Accessories and Staff Education

Developed to assist you in facility protocol compliance, meet regulatory as well as accreditation requirements, and provide a comprehensive training program to ensure proper use. Visit pdihc.com for more education.



PDI 3-in-1 Universal Sani-Bracket®
for wall/pole mounting of a Sani-Cloth® Canister



Sani-Canister Caddy®
for placement on flat surfaces (bilingual)



Sani-TAG™ Equipment ID System
for identification of equipment to ensure the proper wipe is used on a particular piece of equipment




Instructions for Use Signage
(bilingual)



Quick Reference Guide
(bilingual)



Daily Cleaning Wall Chart
(bilingual)



Terminal Cleaning Wall Chart
(bilingual)



High-Touch Surface Guide
(bilingual)



Sani-Cloth® Bleach Germicidal Disposable Wipe

Ideal for disinfecting high-risk areas endemic with *Clostridium difficile* spores and norovirus.

Large Canister	P54072	6 x 10.5"	12/75s	26.25 lbs	1.54 ft	10/3
X-Large Canister	P25784	7.5 x 15"	6/65s	20.50 lbs	1.01 ft	10/4
Large Individual Packet	H58195	5 x 7"	10/40s	8.10 lbs	0.51 ft	19/5
X-Large Individual Packet	U26595	11.5 x 11.75"	3/40s	8.20 lbs	0.50 ft	24/3
Pail	P7007P	7.5 x 15"	2/160s	16.50 lbs	0.95 ft	10/5
Refill (for Pail)	P700RF	7.5 x 15"	2/160s	15.20 lbs	0.53 ft	18/3



Sani-Cloth® AF3 Germicidal Disposable Wipe

Protecting your staff, patients, and residents has never been easier! Alcohol and fragrance-free formula is ideal for use around those with respiratory sensitivities.

Large Canister	P13872	6 x 6.75"	12/160s	26.00 lbs	1.546 ft	10/3
X-Large Canister	P63884	7.5 x 15"	6/65s	15.30 lbs	1.01 ft	10/4
Large Individual Packet	H59200	5 x 8"	10/50s	7.72 lbs	0.54 ft	20/4
X-Large Individual Packet	U27500	11.5 x 11.75"	3/50s	8.00 lbs	0.50 ft	24/3
Portable Pack	M8063S80	8.2 x 9.8"	9/80s	18.05 lbs	0.95 ft	10/4
Pail	P1450P	7.5 x 15"	2/160s	12.40 lbs	0.95 ft	10/5
Refill (for Pail)	P2450P	7.5 x 15"	2/160s	11.10 lbs	0.53 ft	18/4



Super Sani-Cloth® Germicidal Disposable Wipe

Fast, 2-minute contact time allows for a quick room turnover. Ideal for daily use in fast-paced environments that require short contact times and broad coverage of microorganisms.

Large Canister	Q55172	6 x 6.75"	12/160s	26.00 lbs	1.55 ft	10/3
X-Large Canister	Q86984	7.5 x 15"	6/65s	15.30 lbs	1.01 ft	10/4
Large Individual Packet	H04082	5 x 8"	10/50s	7.72 lbs	0.539 ft	20/4
X-Large Individual Packet	U87295	11.5 x 11.75"	3/50s	8.00 lbs	0.50 ft	24/3



Sani-Cloth® Prime Germicidal Disposable Wipe & Sani-Prime® Germicidal Spray

Powered by ACCELOQUAT™ formulation, protect your patients, staff, and facility from HAIs in just 1 minute.

Large Canister	P25372	6 x 6.75"	12/160s	16.00 lbs	1.410 ft	10/3
X-Large Canister	P24284	7.5 x 15"	6/70s	27.84 lbs	1.010 ft	10/4
Large Individual Packet	H06182	5 x 8"	10/50s	7.99 lbs	0.540 ft	20/4
X-Large Individual Packet	U13195	11.5 x 11.75"	3/50s	8.65 lbs	0.500 ft	24/3
Spray Bottle	X12309	n/a	9/32oz per bottle	19.24 lbs	0.828 ft	12/4



Easy Screen® Cleaning Wipe

A compatible touchscreen cleaner that works at the speed of technology. Contains 70% IPA.

Medium Canister	P03672	6" X 9"	12/70s	10.64 lbs	.990 ft	8/5
-----------------	--------	---------	--------	-----------	---------	-----

POINT-OF-CARE ACCESSORIES



3-in-1 Universal Sani-Bracket®

Attach **Sani-Cloth** and **Sani-Hands** canisters to walls or equipment for convenient access; customize with product specific labels.

3-in-1 Universal Sani-Bracket (10 pack)	P010900	n/a	10	19.50 lbs	2.03 ft	4/4
3-in-1 Universal Sani-Bracket (1 pack)	P010901	n/a	1	1.75 lbs	0.17 ft	48/5
Sani-Bracket for X-Large Box	P44210	n/a	1	0.70 lbs	0.07 ft	72/8
3-in-1 Sani-Bracket Pole Accessory Kit	P010600	n/a	1	1.65 lbs	0.05 ft	63/14
Bracket Adapter	P015300	n/a	10	0.23 lbs	n/a	n/a



Sani-Canister Caddy®

Countertop placeholder for canisters. Displays important product and compliance information.

Sani-Canister Caddy Sani-Hands (10 pack)	P013500	n/a	10	1.50 lbs	0.25 ft	28/6
Sani-Canister Caddy Sani-Hands (5 pack)	P013501	n/a	5	0.80 lbs	0.11 ft	30/10
Sani-Canister Caddy Super Sani-Cloth (10 pack)	P012600	n/a	10	1.50 lbs	0.25 ft	28/6
Sani-Canister Caddy Super Sani-Cloth (5 pack)	P012601	n/a	5	0.80 lbs	0.11 ft	30/10
Sani-Canister Caddy Sani-Cloth AF3 (10 pack)	P015200	n/a	10	1.50 lbs	0.25 ft	28/6
Sani-Canister Caddy Sani-Cloth AF3 (5 pack)	P015201	n/a	5	0.80 lbs	0.11 ft	30/10
Sani-Canister Caddy Sani-Cloth Bleach (10 pack)	P013100	n/a	10	1.50 lbs	0.25 ft	28/6
Sani-Canister Caddy Sani-Cloth Bleach (5 pack)	P013101	n/a	5	0.80 lbs	0.11 ft	30/10
Sani-Canister Caddy Easy Screen	P016101	n/a	5	0.80 lbs	0.11 ft	30/10



Sani-Cloth® AF3 Portable Pack Bracket
 Attach the portable pack to mobile equipment for ease of use and accessibility.

Sani-Cloth AF3 Portable Pack Mobile Equipment Bracket	P013200	n/a	1	0.10 lbs	n/a	n/a
Sani-Cloth AF3 Portable Pack Bracket	P039100	n/a	10	n/a	n/a	n/a



Sani-Bracket® Deluxe Multi-Canister Floor Stand
 Place in entryways and other high traffic areas. Holds up to 5 canisters.

Sani-Bracket Deluxe Multi-Canister Floor Stand	P013600	n/a	1	41.00 lbs	8.01 ft	2/2
------------------------------------------------	---------	-----	---	-----------	---------	-----



PDI Safety Stand
 Holds one canister of product and has optional attachments that can be affixed directly to the floor stand to hold personal protective equipment.

PDI Safety Stand	P013700	n/a	1	18.00 lbs	1.81 ft	3/7
PPE Station	P014100	n/a	1	5.50 lbs	0.64 ft	11/7



Sani-TAG® Equipment ID System
 Adhere to equipment to help staff easily identify it as clean. Removes easily and does not leave a sticky residue. Meets regulatory and accreditation standards for identifying clean vs. dirty equipment.

Sani-TAG-Equipment ID System (Clean Tags)	P012200	n/a	8/100s	5.50 lbs	0.29 ft	12/5
-------------------------------------------	---------	-----	--------	----------	---------	------

*Sani-TAG - Disinfect with Sani-Cloth Bleach	*Sani-TAG - Disinfect with Sani-Cloth AF3
*Sani-TAG - Disinfect with Sani-Cloth Plus	*Sani-TAG - Disinfect with Sani-Cloth HB
*Sani-TAG - Disinfect with Super Sani-Cloth	*Sani-TAG - Please use Easy Screen Cleaning Wipe

*Request from sales representative

For more information, or to place an order, contact your local PDI sales representative.

#1 DISINFECTANT WIPE IN HEALTHCARE¹



¹GHX 2018 Q3 report.